



PLANCHETS and RADIO CHEMISTRY DISCS



Type	P/N		Diameter Overall		Height Overall		Recessed Dia.	
	Aluminum	Stainless Steel	Inches	MM	Inches	MM	Inches	MM
Cupped	C-1A	C-1S	1	25.4	3/32	2.4		
	C-2A	C-2S	1	25.4	1/4	6.4		
	C-3A	C-3S	1	25.4	5/16	7.9		
	C-4A	C-4S	1-1/4	31.8	3/32	2.4		
	C-5A	C-5S	2	50.8	3/32	2.4		
	C-5AA	C-5AS	2	50.8	1/8	3.2		
	C-6A	C-6S	2	50.8	5/32	4.0		
	C-7A	C-7S	2-1/16	52.4	3/32	2.4		
	C-8A	C-8S	2-1/16	52.4	5/32	4.0		
	C-9A	C-9S	2	50.8	3/16	4.8		
	C-10A	C-10S	2	50.8	1/4	6.4		
	C-11A	C-11S	2	50.8	5/16	7.9		
	C-12A	C-12S	2-3/8	60.3	1/8	3.2		
	C-12AA	C-12AS	3-1/8	80.0	7/32	5.5		
	C-12BA	C-12BS	2-3/8	60.3	1/4	6.4		
C-13A	C-13S	4-3/16	106.4	1/4	6.4			
C-13AA	C-13AS	4-3/16	106.4	1/8	3.2			
C-14A	C-14S	4-1/2	114.3	1/8	3.2			
C-15A	C-15S	4-7/8	123.8	5/32	4.0			
Dished	D-AA	D-AS	1	25.4	3/32	2.4		
	D-1A	D-1S	2	50.8	5/32	4.0		
	D-2A	D-2S	2-1/16	52.4	5/32	4.0		
Flanged	FL-1A	FL-1S	1-5/8	41.3	5/64	2.0	1-5/16	33.3
	FL-2A	FL-2S	1-5/8	41.3	1/8	3.2	1	25.4
	FL-3A	FL-3S	2-5/16	58.9	5/32	4.0	2-3/32	53.2
	FL-4A	FL-4S	2-5/16	58.9	5/16	7.9	2-3/32	53.2
Flat	F-CA	F-CS	*A	5/8	15.9	*PACKAGING PREFERENCE DESIGNATION: A, C OR D: ALL POLISHED ONE SIDE A-BULK C-WITH PROTECTIVE TAPE D-ONE PER ENVELOPE		
	F-BA	F-BS		11/16	17.5			
	F-AA	F-AS	*C	3/4	19.0			
	F-AAA	F-AAS		7/8	22.0			
	F-1A	F-1S		*D	1		25.4	
	F-2A	F-2S	1-1/4		31.8			
	F-3A	F-3S	1-1/2		38.1			
	F-3AA	F-3AS	1-3/4		44.0			
F-4A	F-4S	2	50.8					

SAMPLES ON REQUEST

MAY 2019


PLANCHETS and RADIOCHEMISTRY DISCS

SPECIAL ORDER


TYPE	STOCK NUMBER	DIAMETER OVERALL	
	NICKEL*	INCHES	MM
FLAT	F-CN	5/8	15.9
	F-BN	11/16	17.5
	F-AN	3/4	19.0
	F-AAN	7/8	22.0
	F-IN	1	25.4

*Other precious metals : Silver, Platinum, Copper, Gold, Tantalum, etc., also available

PLANCHETS & RADIOCHEMISTRY DISCS (CONCENTRIC BOTTOM)



TYPE	STOCK # ALUMINUM	STOCK # STAINLESS STEEL	DIAMETER OVERALL		HEIGHT OVERALL	
			INCHES	MM	INCHES	MM
CUPPED	C-1AC	C-1SC	1	25.4	3/32	2.4
	C-2AC	C-2SC	1	25.4	1/4	6.4
	C-3AC	C-3SC	1	25.4	5/16	7.9
	C-4AC	C-4SC	1-1/4	31.8	3/32	2.4
	C-5AC	C-5SC	2	50.8	3/32	2.47
	C-5AAC	C-5ASC	2	50.8	1/8	3.2
	C-6AC	C-6SC	2	50.8	5/32	4.0
	C-7AC	C-7SC	2-1/16	52.4	3/32	2.4
	C-8AC	C-8SC	2-1/16	52.4	5/32	4.0
	C-9AC	C-9SC	2	50.8	3/16	4.8
	C-10AC	C-10SC	2	50.8	1/4	6.4
	C-11AC	C-11SC	2	50.8	5/16	7.9
	C-14AC	C-14SC	4-1/2	114.3	1/8	3.2



MARINELLI BEAKER SELECTION CHART BASED UPON DETECTOR ENDCAP DIAMETER

Detector Endcap Diameter =	2.30" (5.8cm)	2.50" (6.4cm)	2.75" (7.0cm)	3.00" (7.6cm)	3.15" (8.0cm)	
Available Models :	523N-E 523GWMTJ	125G-E G-125G	527G-E 127G-E 227G-E G-127G 127GWMTJ	443016 530G-E 130G-E 230G-E 430G-E G-530G G-130G G-430G 130GWMTJ 530GWMTJ 530GSWMTJ	580G-E	
Detector Endcap Diameter =	3.15"(8.0cm) or 3.25"(8.3cm)	3.54" (9.0cm)	3.54"(9.0cm) or 3.75"(9.5cm)	4.00" (10.2cm)	4.25" (10.8cm)	4.75" (12.1cm)
Available Models :	463316 533N-E 132G-E 133N-E 233N-E 433N-E G-133N G-433N 133GWMTJ	590G-E 190G-E	538G-E 138G-E G-438G	541G-E 141G-E 441G-E G-441G	445N-E G-445N	448G-E

Fast & Friendly Customer Service

Ordering by phone
(352) 687-8840



Ordering by Email
gama@ga-maassociates.com

Ordering by Fax
(352) 687-2117



For your convenience
MasterCard, Visa, American Express & Wire Transfer

For Technical Assistance
(352) (687-8840)

